REF 022 500 SP

# EpiContur Sport

**WORK & PLAY**

Made in  
Germany

## SPRACHEN LANGUAGES

<b>DE</b>	Deutsch	Gebrauchsanweisung .....	04
<b>EN</b>	English	Instructions for use .....	08

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## EN – EpiContur Sport

Many thanks for placing your trust in this medical device. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### INTENDED PURPOSE

This medical device is a brace for relief of the muscle origins on the ulnar and/or radial epicondyle.

### INDICATIONS

Epicondylitis, chronic or recurring epicondylopathy, tendopathies, excessive strain, so-called tennis/golf arm.

### CONTRAINDICATIONS

Peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders with the region of the body treated, skin diseases in the part of the body treated.

### APPLICATION RISKS / IMPORTANT NOTES

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e. g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds

- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

### **FITTING INSTRUCTIONS**

By loosening the Velcro fastening in the fitted central area, both parts of the brace can be pushed together. Therefore, the brace can be set to the individual circumference from 21 cm to 35 cm steplessly. Here, ensure that the two widened areas of the brace attached on the upper and lower side of the arm are placed as close as possible to opposite each other. The pad can be placed freely using the Velcro fastener. For correct positioning, observe the markings printed on the support. The narrower area of the pad should point towards the elbow joint. Thread the Velcro tape through the plastic eyelet and then position the brace on the lower arm. The pad should exert pressure on the muscle belly which leads from the hand to the elbow. The

muscle belly can be clearly felt when you make a fist. The EpiContur brace should only exert pressure on this muscle if you make a fist or grasping movement. If the hand is relaxed, the brace should be placed on the pad without noticeable pressure. Close the Velcro fastener so that the brace fits closely but causes no congestion.

### **REMOVAL**

To remove the epicondylitis brace, loosen the Velcro tape, increase the circumference and remove the brace in the direction of the hand.

### **MATERIAL COMPOSITION**

Polyamide (PA), polyurethane (PU), polyester (PES)  
The sewn in textile label provides the precise material composition.

## CLEANING INFORMATION

 Hand-washing 30°C  Do not bleach

 Do not dry in a tumble dryer  Do not iron

 Do not clean chemically

Do not use fabric conditioner. Dry in the fresh air.

## GUARANTEE

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

## USEFUL LIFE / LIFETIME OF THE PRODUCT

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

## DUTY OF NOTIFICATION

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find the manufacturer's contact data in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

## DISPOSAL

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

## DECLARATION OF CONFORMITY

The manufacturer confirms that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL.

Status: 04.2024

**MD** Medical device |  Single patient – multiple use

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**VPS Medical Ltd.**
**BORT GmbH**